
Experiencing & Utilizing Value-Adding Features of FINO.VITA GX PACS for Higher Vocational Education on Diagnostic Imaging

Picture Archiving and Communication System (PACS) is a system to store, archive, view and distribute images across hospital departments.

PACS has been well known and used in the past few decades globally, however in Indonesia, this has just become a trend due to the introduction of full Electronic Medical Record and SATU SEHAT Health Technology Platform. This platform is one of the 6 pillars¹ of Indonesia Health Transformation that are prepared to assist on managing the 9 diseases priorities².

Transformasi sistem kesehatan 2021-2024

5 RPJMN dan 6 pilar transformasi

The platform will integrate all patient health data across all hospitals so that it could be accessible for physicians and patients alike. This initiative is meant to make access to patient data easier to enable holistic and informed decision making.

As a trend, this invites many PACS and Healthcare IT companies both international and local corporations to become relevant players here in Indonesia.

On this brief testimonials, we will address the importance of Medical Device Certified DICOM Viewer of a PACS to help hospitals select the right PACS for their doctors and patients and how FINO.VITA GX could give higher value-add to them in response to managing the 9 diseases priorities above.

¹Primary Care, Referral Service, Health Resilience, Healthcare Financing, Healthcare Human Resources, Health Technology

²Cancer, Cardiovascular, Stroke, Gastrohepatology, Uronefrology, Mother and Child Health, Tuberculosis, Diabetes Melitus, and Infection Diseases.

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
As a Radiologist Viewer, DICOM Viewer usually consists of some tools that include basic measurement as well as image processing tools. These tools have been tested for accuracy and quality assurance (QA) to ensure consistent image interpretation.

However, **many free and/or DICOM Viewers we found easily on internet are not through this testing and QA and even have put a disclaimer that the Viewer is not to be used for Diagnostic Purposes.** Thus, this can lead to misdiagnosis and missed diagnosis, that will put multiple stakeholders at risk, especially Patients and Reading Physicians.

As a commercial PACS, Konica Minolta provides more value-add to hospitals, physicians and patients with following features within the Viewer :

1. Bone Suppression
2. Integrated Reporting Workflow (FINO.Report)

BONE SUPPRESSION

“By comparing Chest x-ray using Konica Minolta *Bone suppression images* with Chest x-ray standard on patients with TB MDR. The overall result is showing that the *Bone Suppression Images* is effective to helping radiologist analyse all abnormalities in the entire lung area, including the apex area. *Bone Suppression Images* eliminate images of rib/costae so chest area can be evaluated clearly without superposition between rib/costae with the soft tissue of the lung. *Bone suppression images* from Konica Minolta also can provides more information from chest x-ray and doctor can resulting the accurate diagnosis.”³


In addition, the feature could also help physicians to spot nodule better that could later be combined with AI – Chest Xray Nodule and Tuberculosis Detection.

¹Primary Care, Referral Service, Health Resilience, Healthcare Financing, Healthcare Human Resources, Health Technology

²Cancer, Cardiovascular, Stroke, Gastrohepatology, Urology, Mother and Child Health, Tuberculosis, Diabetes Melitus, and Infection Diseases.

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INTEGRATED REPORTING WORKFLOW (FINO.Report)

A RIS-driven and PACS-driven workflow are 2 workflows that are available when it comes to reporting findings and diagnosis of a radiology case. PACS-driven workflow allows user to open the study first followed by opening the reporting/expertise. This allows users to avoid report-study mismatch that might happen in RIS-driven workflow. In addition, an integrated reporting flow allows a single-database running for both studies, reports and any other information, making accessibility to any patient information a lot easier.

Secondly, due the ability to compare past and current report and then seamlessly cite previous report, this would allow better, seamless and holistic decision-making processes for radiologists, hence ensuring better management of diseases, per SATUSEHAT intention. Finally, the ability to attach key images on the expertise body will help boost clarity for Referring Physicians.

In conclusion, both features as described above could potentially support both Digital Transformation and Health Transformation of Indonesia, emphasizing on assisting the detection of lung nodule and tuberculosis as well as supporting seamless and safe reporting workflow.

Dr. Gatot Murti Wibowo, S.Pd., M.Sc.

Ketua Program Studi Magister Imaging Diagnostik Pascasarjana

Poltekkes Kemenkes Semarang

¹Primary Care, Referral Service, Health Resilience, Healthcare Financing, Healthcare Human Resources, Health Technology

²Cancer, Cardiovascular, Stroke, Gastrohepatology, Urology, Mother and Child Health, Tuberculosis, Diabetes Mellitus, and Infection Diseases.

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